




The **CONTACT** family is complete.



CONTACTLINSEN

 Biohydrogel

 Gentle Edge Profile



 Made in Germany

A brand of  wöhlk CONTACTLINSEN



LIFE.

SPHERIC

CONTACT LIFE. SPHERIC

1. Spherical monthly lenses with smooth edge profile
2. Comfortable for all-day wear
3. Great option for patients who experience dry eye or discomfort with contact lenses
4. Designed for lipid-rich lacrimal fluid
5. Ideal for screen-based work

TORIC

CONTACT LIFE. TORIC

1. Toric monthly lenses with smooth edge profile
2. Comfortable for all-day wear
3. Great option for patients who experience dry eye or discomfort with contact lenses
4. Designed for lipid-rich lacrimal fluid
5. Ideal for screen-based work
6. Provides stable vision for astigmatism
7. Available in -2.75D cylinder across all axes

MULTI SPHERIC


CONTACT LIFE. MULTI SPHERIC

1. Distance-centred multifocal spherical monthly lenses with smooth edge profile
2. Comfortable for all-day wear
3. Great option for patients who experience dry eye or discomfort with contact lenses
4. Designed for lipid-rich lacrimal fluid
5. Ideal for screen-based work
6. Corrects presbyopia

MULTI TORIC

CONTACT LIFE. MULTI TORIC


1. Distance-centred multifocal toric monthly lenses with smooth edge profile
2. Comfortable for all-day wear
3. Great option for patients who experience dry eye or discomfort with contact lenses
4. Designed for lipid-rich lacrimal fluid
5. Ideal for screen-based work
6. Stable vision for presbyopia with astigmatism



SPHERIC

Availability CONTACT LIFE. SPHERIC

BASE CURVE [mm]	DIAMETER [mm]	SPHERE [D]
8.30	13.60	-12.00 to -8.00 (0.50) -8.00 to +5.00 (0.25) +5.00 to +8.00 (0.50)
8.60		-20.00 to -8.00 (0.50) -8.00 to +5.00 (0.25) +5.00 to +10.00 (0.50)
8.80 (+)	14.20	plan to +5.00 (0.25) +5.00 to +8.00 (0.50)
8.90 (-)		-12.00 to -8.00 (0.50) -8.00 to plan (0.25)




steps 0.25 steps 0.50

Price

Box of 6	€ 21.78
Box of 3	€ 13.00
Blister *2	€ 2.50


*2 Maximum order quantity = 2 blisters with the same parameters.



TORIC

Availability CONTACT LIFE. TORIC

BASE CURVE [mm]	DIAMETER [mm]	SPHERE [D]
8.60	14.40	-8.00 to -6.00 (0.50) -6.00 to +4.00 (0.25) +4.00 to +6.00 (0.50)




steps 0.25 steps 0.50

CYLINDER [D]	-0.75 to -2.75 (0.50)
AXIS [°]	10 to 180 (10)

Price

Box of 6	€ 39.05
Box of 3	€ 22.00
Blister *2	€ 2.50

*2 Maximum order quantity = 2 blisters with the same parameters.



MULTI SPHERIC

N = Near / D = Distance

Availability CONTACT LIFE. MULTI SPHERIC

BASE CURVE [mm]	DIAMETER [mm]	SPHERE [D]
8.60	14.20	-8.00 to +6.00 (0.25)

steps 0.25


ADDITION [D]	Typ A (+1.00) / Typ B (+1.75)	(0.75)

Price

Box of 6	€ 48.50
Box of 3	available soon
Blister *2	€ 2.50

*2 Maximum order quantity = 2 blisters with the same parameters.

Available soon



MULTI TORIC

N = Near / D = Distance

Availability CONTACT LIFE. MULTI TORIC

BASE CURVE [mm]	DIAMETER [mm]	SPHERE [D]
8.60	14.40	-5.00 to +3.00 (0.25)

steps 0.25

CYLINDER [D]	-0.75 / -1.25

AXIS [°]	20/70/90/110/160/180

ADDITION [D]	Typ A (+1.00) / Typ B (+1.75)	(0.75)

Price

Box of 6	available soon
Box of 3	available soon
Blister *2	available soon

*2 Maximum order quantity = 2 blisters with the same parameters.




Dear contact lens specialists,

It's finally here: the multifocal version of our popular **CONTACT LIFE** monthly lens. The spherical and toric versions are now being extended by the '**MULTI SPHERIC**' version and will also be available as '**MULTI TORIC**' in a few months' time – making our **CONTACT LIFE** family complete!

The gentle edge profile has been carried forward into the new LIFE versions as well. It is slightly rounded and fulfils two functions at once: The lacrimal fluid under the contact lens is optimally exchanged and the eye is better supplied with oxygen.

Edge Profile



- ▶ Minimal influence on the physiological processes between the conjunctiva and bulbar conjunctiva
- ▶ Optimum exchange of lacrimal fluid
- ▶ Long-term compatibility guaranteed

- ▶ Creation of edge imprints and scuff marks
- ▶ Reduced exchange of lacrimal fluid
- ▶ Reduced long-term compatibility

Material

Our monthly multifocal lenses are also made from the popular and proven **Vitafilcon A** bio-hydrogel. It is excellent at retaining moisture and is highly resistant to deposits. The Vitafilcon material was developed by us and was honoured with the Sicca Award by the German Ophthalmologists' Association in 2004!

Our Material has been honoured with an award!

with the



SICCA-FÖRDERPREIS

Made in Germany

As both the Vitafilcon material and the monthly disposable lenses are manufactured in Schönkirchen, our **CONTACT LIFE** is 100 % made in Germany. We hope you enjoy testing our new monthly multifocal lenses and, above all, that you have many successful fittings. If you have any questions, please contact us at any time!



CONTACTLINSEN

WÖHLK CONTACTLINSEN GMBH

ADDRESS:

Bürgermeister-Schade-Str. 12-16
24232 SCHÖNKIRCHEN / GERMANY

SERVICE HOURS:

Mon - Fri 9:30 a.m - 3:00 p.m.

SERVICE EXPORT:

Phone +49(0)431 991-12 00
Fax +49(0)431 991-15 55
Email export@wohlc.com
Internet: wohlc.com

Art.7826 /REV.01 /Stand 14.01.2026



A brand of  **wohlc** CONTACTLINSEN