

# PRODUCT OVERVIEW






**wöhlk**  
CONTACTLINSEN



# ZEISS Contact Day 1

Daily contact lenses



PRODUCT	DESCRIPTION	MATERIAL	GEOMETRY	AVAILABLE PARAMETERS [mm]	POWERS [D]
<b>ZEISS Contact Day 1</b> SPHERIC <a href="#">Art.Nr. 2910 (Box of 32)</a> 	Spheric <b>ZEISS</b> daily replacement lens from bio-hydrogel for correcting spheric ametropias	Bio-hydrogel 58% water bi-ionic Tint: light blue Dk/t (-3.0): $42.9 \times 10^{-11}$ Barrer	<b>Back surface:</b> bi-curved spheric with smooth edge profile <b>Front surface:</b> spheric front optical zone (VOZ) with lenticular transition zone, adaptive design control	BC: 8.8 DIA:14.20	sphere: -16.00 to -6.50 (0.50) -6.00 to -0.50 (0.25) plan +0.50 to +5.00 (0.25) +5.50 to +8.00 (0.50)
<b>ZEISS Contact Day 1</b> TORIC <a href="#">Art.Nr. 2941 (Box of 32)</a> 	Bi-curved toric <b>ZEISS</b> daily replacement lens with smooth edge profile for correcting astigmatic ametropias	Bio-hydrogel 58% water bi-ionic Tint: light blue Dk/t (-3.0): $27.3 \times 10^{-11}$ Barrer	<b>Back surface:</b> bi-curved toric with smooth edge profile <b>Front surface:</b> dynamic prism ballast, upper slab-off, adaptive design control	BC: 8.8 DIA:14.20	sphere: -10.00 to -6.50 (0.50) -6.00 to plan (0.25) cylinder: -1.75 / -1.25 / -0.75 (-2.25 in 180°) axis [°]: 10/20/70/90/110/160/180 <b>Parametererweiterung</b> sphere: +0.50 to +5.00 (0.50) cylinder: -1.75 / -1.25 / -0.75 axis [°]: 20/90/160/180
<b>ZEISS Contact Day 1</b> MULTIFOCAL <a href="#">Art.Nr. 2923 (Box of 32)</a> 	Multifocal <b>ZEISS</b> daily wear lens for correcting spheric age-related ametropias	Bio-hydrogel 58% water bi-ionic Tint: light blue Dk/t (-3.0): $42.9 \times 10^{-11}$ Barrer	<b>Back surface:</b> bi-curved spheric with smooth edge profile <b>Front surface:</b> multifocal aspheric distance-centered front optical zone (VOZ) with lenticular transition zone, adaptive design control	BC: 8.8 DIA:14.20	sphere: -10.00 to +5.00 (0.25) addition: „Type A“ (+0.75) „Type B“ (+1.50)



# ZEISS Contact Day 30

Monthly contact lens



NEW

Coming soon!



wöhlk  
CONTACTLINSEN

PRODUCT	DESCRIPTION	MATERIAL	GEOMETRY	AVAILABLE PARAMETERS [mm]	POWERS [D]
<b>ZEISS Contact Day 30 compatic SPHERIC</b> Art.Nr. 2640 (Box of 6)	Monthly replacement lens made of bio-hydrogel with smooth edge profile for comfortable daily wear.	Vitafilcon A Tint: light green 53% water bi-ionic neutral	bi-curved spheric	BC: 8.3 DIA: 13.60	sphere: -12.00 to -8.00 (0.50) -8.00 to +5.00 (0.25) +5.00 to +8.00 (0.50)
				BC: 8.6 (+/-) DIA: 14.20	sphere: -20.00 to -8.00 (0.50) -8.00 to +5.00 (0.25) +5.00 to +10.00 (0.50)
				BC: 8.8 (+) DIA: 14.20	sphere: plan to +5.00 (0.25) +5.00 to +8.00 (0.50)
				BC: 8.9 (-) DIA: 14.20	sphere: -12.00 to -8.00 (0.50) -8.00 to plan (0.25)
<b>ZEISS Contact Day 30 compatic TORIC</b> Art.Nr. 2401 (Box of 6)	Toric monthly replacement lens made of bio-hydrogel with smooth edge profile and dynamic symmetric stabiliser for comfortable daily wear.	Vitafilcon A Tint: light green 53% water bi-ionic neutral	bi-curved back surface toric dynamic symmetric	BC: 8.6 DIA: 14.40	sphere: -8.00 to +6.00 (0.50) -6.00 to +4.00 (0.25) cylinder: -0.75 to -2.75 (0.50) axis [°]: 10 to 180 (10)
<b>ZEISS Contact Day 30 compatic MULTIFOCAL SPHERIC</b> Art.Nr. 2643 (Box of 6)	Monthly distance-centered multi-focal spherical lens made of bio-hydrogel with gentle, double-sided edge profile for comfortable daily wear.	Vitafilcon A Tint: light green 53% water bi-ionic neutral	bi-curved spheric	BC: 8.6 DIA: 14.20	sphere: -8.00 to +6.00 (0.25) addition: Typ A (+1.00) Typ B (+1.75) (0.75)
<b>ZEISS Contact Day 30 compatic MULTIFOCAL TORIC</b> Art.Nr. 2403 (Box of 6)	Monthly distance-centered multifocal torical lens made of bio-hydrogel with gentle, double-sided edge profile for comfortable daily wear.	Vitafilcon A Tint: light green 53% water bi-ionic neutral	bi-curved back surface toric dynamic symmetric	BC: 8.6 DIA: 14.40	sphere: -5.00 to +3.00 (0.25) cylinder: -0.75 / -1.25 (0.50) addition: Typ A (+1.00) Typ B (+1.75) (0.75) axis [°]: 20/70/90/110/160/180



# ZEISS Care Solution



PRODUCT	DESCRIPTION	CONTENTS	ART-NO.	UNIT	AMOUNTE
<b>ZEISS All-in-one ADVANCE</b>	Extremely compatible and effective all-in-one <b>ZEISS</b> solution for all soft contact lenses.	Solution CL container Supplementary sheet	7302	24	360 ml
			7342	12	100 ml



**wöhlk**  
CONTACTLINSEN



**Distribution: wöhlk** Contactlinsen GmbH

**Address:**  
Bürgermeister-Schade-Str. 12-16  
DE - 24232 Schönkirchen

[woehk.com](http://woehk.com)

**Service hours:**  
Mon - Fri 9:30 a.m. - 3:00 p.m.

**Service online:**  
Email [info@woehk.com](mailto:info@woehk.com)  
Internet [woehk.com](http://woehk.com)

**Service Export:**  
Phone +49 (0) 431 991 - 1200  
Fax +49 (0) 431 991 - 1515  
Email [export@woehk.com](mailto:export@woehk.com)

**Distribution Export:**  
Phone +49 (0) 151 510 94 738  
Fax +49 (0) 431 991 -1515  
Email [export@woehk.com](mailto:export@woehk.com)

ZEISS is a registered trademark of Carl Zeiss AG used under license of Carl Zeiss Vision International GmbH.